

JUMPING SPIDER

— Powered by Geyi Medical —

Portable Image Processor

One Tablet, Multiple Endoscopy Solutions.



NMPA



PRODUCT FEATURE

Ureteroscope
Single-use



Cystoscope
Single-use



Choledochoscope
Single-use



Bronchoscope
Single-use



Hysteroscope
Single-use



12.1" Touch Screen

Intuitive touch operation, 1280 × 800 Display / 1920 × 1080 Extended Output.

Compatible With Various Endoscopes

Type-C Connector

Up to 8 Hours Battery

Portable and flexible operation.

32 GB Built-in Storage

Convenient image & video export



Special Features



Photograph



Record



Playback



LED Adjustment



Freeze



White Balance



Image Ratio



Image Zoom in



Sharpness Adjustment



RGB Color Adjustment



Eliminate Image Reflections



Multi-Language Support

PRODUCT PARAMETER

Model	Video Output	Display Resolution	Built-in Battery
GY-AIO-121	HDMI1.4 1920×1080 (1080P)	1280×800 (720P)	Excluded
GY-AIO-121(A)			Included

ZHEJIANG GEYI MEDICAL INSTRUMENT CO., LTD.

✉ info@geyi-medical.com

🌐 www.geyi-medical.com

🌐 www.jumpingspider-med.com