

JUMPING SPIDER

— Powered by Geyi Medical —

Single-use Bronchoscope

Airway, Always Ready



EU-MDR



NMPA



PRODUCT FEATURE

Handle Lever

Uniform damping on handle lever, instant stop at any position.

Type-C Connector

Suction Button

Suction Connector

Custom Button

Ergonomic Handle

With anti-slip surface treatment.

Working Channel

Direct-access instrument channel enables convenient operation.

Rotation Part

90° L & R rotatable

High Sensitivity CMOS+LED Lights


HD Image

Megapixel clear field of view, assists clinical discrimination and observation of bronchial airway tissues.



220°

Bidirectional Bending



Flexible bronchoscope with precise deflection; multiple specifications to meet diagnosis & treatment demands of different patient groups.

PRODUCT PARAMETER

Model	Outer Diameter	Working Channel	Working Length (mm)	Bending Angle	Depth of Field (mm)	Viewing Angle	Resolution (pixel)
GY-BRM075A1L	7.5 Fr	3.3 Fr	600	220° (up & down)	3-50	120°	160,000
GY-BRM084A1L	8.4 Fr	3.6 Fr	600				
GY-BRM120A1L	12.0 Fr	6.0 Fr	600				
GY-BRN147A1L	14.7 Fr	6.6 Fr	600		3-100	140°	1,000,000

ZHEJIANG GEYI MEDICAL INSTRUMENT CO., LTD.

✉ info@geyi-medical.com

🌐 www.geyi-medical.com

🌐 www.jumpingspider-med.com