


JUMPING SPIDER

— Powered by Geyi Medical —

Single-use Choledochoscope


Into the Duct, Into View



NMPA



PRODUCT FEATURE




HD Image

Megapixel clear field of view, assists clinical observation and differentiation of biliary tract tissues.



220°

Bidirectional Bending



Flexible choledochoscope with precise deflection; multiple specifications to meet diagnosis & treatment demands of different patient groups.

PRODUCT PARAMETER

Model	Outer Diameter	Working Channel	Working Length (mm)	Bending Angle	Depth of Field (mm)	Viewing Angle	Resolution (pixel)
GY-CTM075A1BL	7.5 Fr	3.3 Fr	600	220° (up & down)	3-50	120°	160,000
GY-CTM084A1BL	8.4 Fr	3.6 Fr	600				
GY-CTM120A1BL	12.0 Fr	6.0 Fr	600				
GY-CTM120A1BH	12.0 Fr	6.0 Fr	400				
GY-CTN147A1BL	14.7 Fr	6.6 Fr	600		3-100	140°	1,000,000
GY-CTN147A1BH	14.7 Fr	6.6 Fr	400				

ZHEJIANG GEYI MEDICAL INSTRUMENT CO., LTD.

✉ info@geyi-medical.com

🌐 www.geyi-medical.com

🌐 www.jumpingspider-med.com