

JUMPING SPIDER

— Powered by Geyi Medical —

Single-use Cystoscope

One Scope
Two Roles



EU-MDR



NMPA



PRODUCT FEATURE



220° Bidirectional Deflection

U220° / D220° deflection supports controlled navigation and systematic inspection of the bladder.

As Slim as 8.4 Fr

Designed for smooth transurethral access with a slim insertion profile.



120° / 140° Field of View

Wide viewing angle supports efficient inspection of the bladder.

Suction Connector

One Scope, Two Roles

8.4 Fr, diagnostic purpose;
15.6 Fr, therapeutic purpose.

Designed for Comfortable Handling

Ergonomic handle design supports intuitive and responsive scope manipulation.

Large Working Channel

Up to 6 Fr, provides room for compatible urological instruments during interventional procedures.

Rotation Part

90° L & R rotatable

PRODUCT PARAMETER

Model	Outer Diameter	Working Channel	Working Length (mm)	Bending Angle	Depth of Field (mm)	Viewing Angle	Resolution (pixel)
GY-CYM084WJQB145 aka:GY-CYM084	8.4 Fr	3.6 Fr	400	220° (up & down)	3-50	120°	160,000
GY-CYN156WJQA134 aka:GY-CYN156	15.6 Fr	6.0 Fr			3-100	140°	1,000,000

ZHEJIANG GEYI MEDICAL INSTRUMENT CO., LTD.

✉ info@geyi-medical.com

🌐 www.geyi-medical.com

🌐 www.jumpingspider-med.com