

JUMPING SPIDER

— Powered by Geyi Medical —

Single-use Ureteroscope

As Slim as
6.3Fr




NMPA



PRODUCT FEATURE

Slim 6.3 Fr Distal


To support smooth ureteral access and efficient navigation.



Y adapter included
packaged with each scope.

275° Responsive Bidirectional Deflection

Designed to provide the flexibility needed to access challenging anatomy with greater control.



Handle Lever

- 1:1 coaxial steering, more smooth operation.
- Two shortcut keys, support user customization, can be preset with up to 15 functions.

240 g Ergonomic Lightweight Handle


Designed for comfortable grip and intuitive scope manipulation.

3.6 Fr Working Channel

Supports compatible accessories while helping maintain irrigation performance.

Clear Digital Visualization

Digital imaging supports reliable observation throughout the procedure.



PRODUCT PARAMETER

Model	Outer Diameter	Working Channel	Working Length (mm)	Bending Angle	Depth of Field (mm)	Viewing Angle	Resolution (pixel)
GY-UR6.3	6.3 Fr	3.0 Fr	670	275° (up & down)	3-50	120°	160, 000
GY-UR7.5	7.5 Fr	3.4 Fr					
GY-UR8.4	8.4 Fr	3.6 Fr					
GY-UR9.3	9.3 Fr	3.6 Fr					

ZHEJIANG GEYI MEDICAL INSTRUMENT CO., LTD.

✉ info@geyi-medical.com

🌐 www.geyi-medical.com

🌐 www.jumpingspider-med.com